## Instructions For Use [EN]

Acc. to IEC 82304-1 Health software – Part 1: General requirements for product safety

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### Product Description

Medical product OnyxCeph<sup>3</sup>™ represents a client/server Windows™ application with SQL database for working with 2D and 3D image datasets.

OnyxCeph<sup>3</sup>™ was developed as an image management and image evaluation software for dental units and primarily supports digital workflows of image-based evaluation, treatment planning and patient education in general dentistry, orthodontics and orthognathic oral surgery.

++D981OnyxCeph3S

Release 3.2.211 (585)

### Medical Purpose | Indication for Use

OnyxCeph<sup>3</sup>™ software is intended to be used for the medical purpose of managing and evaluating two-dimensional and three-dimensional images in the context of dental applications by qualified staff only. Clinical decisions cannot be motivated exclusively or even mainly on evaluation results provided by the software.

The software can be used to digitally perform certain image-based orthodontic workflows in addition to metric and angular evaluation of image data.

CH REP

Importer | CH-REP Switzerland

Jordi Röntgentechnik AG  
Eichweid 8a  
6203 Sempach Station  
Schweizerland  
Tel. +41 61 417 9393  
E-Mail: roentgen@rxtech.ch

### R<sub>x</sub> Only

In the U.S., OnyxCeph<sup>3</sup>™ is classified as a Class 2 medical device based upon United States Food and Drug Administration, classification in Orthodontic Software.

Caution: Federal law restricts this device to sale by or on the order of a physician. The product is Rx Only labeled.

## Warnings

- The user must have sufficient experience in using and handling computers as well as Windows™-based software including the Windows™ operating system used.
- The user may only use the software if he/she is proficient in the language selected for the OnyxCeph<sup>3</sup>™ client.
- Use and configure appropriate security measures such as anti-virus software and firewalls to protect against data loss due to hacker attacks.
- Use non-activated trial versions for learning purposes only but not for clinical cases.
- It is the user's responsibility to install available updates in order to always use the latest version.
- Only use OnyxCeph<sup>3</sup>™ if the 2D and 3D image information to be used for the intended workflow can be clearly and verifiably assigned to the correct patient and has the required quality and resolution.

## Notes

- OnyxCeph<sup>3</sup>™ may only be used by qualified personnel like dentists, dental assistants or dental technicians with professional experience in orthodontics.
- Use OnyxCeph<sup>3</sup>™ for clinical cases only if you have been sufficiently instructed in the use of the software for the intended workflow.
- Note that the use of OnyxCeph<sup>3</sup>™ software for digital orthodontic workflows neither influences nor contradicts conventional clinical treatment concepts and treatment principles.
- Do not give your user access data to third parties.
- Please read the manual carefully.
- Any serious incident occurring in connection with the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established (in Germany BfArM).

## Contraindications

There are no known contraindications to the OnyxCeph<sup>3</sup>™ standalone medical software.

## Performance Data

The correct function of the software requires the use of the most up to date program version (release). The availability of a new program version is displayed by a release update button top right in the program main window if online connected. Alternatively, the update information has to be checked on the product website [www.onyxceph.com](http://www.onyxceph.com) on a regular base.

## Accuracy and Clinical Benefit

The functions provided by the software according to its medical purpose are based on the measurement of image data sets. All measurements supported by the application are based on digitized landmark coordinates within the coordinate system of the analyzed image data set while a correct calibration of the image data itself is required.

The accuracy of the measurement of distances and angles on 2D and 3D image data is a minimum of 1/100 of the respective measured quantity.

## Installation

The installation of the software OnyxCeph<sup>3</sup>™ is performed by the user or by its IT partners. To install the software correctly, the setup instructions enclosed in each installation package must carefully followed. For installation, administrator privileges are required and an up to date Windows™ operating system has to be used. Detailed information on installation and configuration is available at [www.onyxwiki.net](http://www.onyxwiki.net). The intended medical purpose requires a correct registration and activation of the software by the user. Unauthorized TRIAL versions or versions labeled TEST can only be used to become familiar with the application acc. to its user interface and functionality but do not represent a medical product.

## Port Sharing

To allow communication between server database and client application, IP port 16300 has to be open on both ends inbound and outbound.

## Start, Application, Exit

After proper installation of the latest program version the main program can be launched on a client computer from the program group OnyxCeph3 or from the related desktop link by double-clicking. Notes and instructions for using the software are provided via the main menu entry [Help] and the context-sensitive help. The program is terminated via the main menu entry [Exit].

## Uninstallation

The workstations are uninstalled using the uninstallation program in the OnyxCeph3 program group. The database server is uninstalled via the Windows <sup>™</sup> control panel in the [Programs and Features] category.

17-11-2023

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gemäß IEC 82304-1 Health software – Part 1: General requirements for product safety

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