
Required Information for Creating a Customized OnyxCeph³™ Analysis

The following information is required to create a customized analysis:

0. Tracing type of the analysis (Lateral Ceph, Pano, Profile Photo, Frontal Foto,)
1. Name of the analysis
2. List of all numeric results to be calculated, incl.
 - 2.1 Short name of the result variable
 - 2.2 Long name of the result variable
 - 2.3 Definition for the result variable incl. no. of digits after decimal point for the numerical value
 - 2.4 Clinical norm (if available) *
 - 2.5 Standard deviation for the norm (if available) *
 - 2.6 Verbal result description for result interbvals (if available)
3. A sketch for the graphical tracing
4. A sketch for the print form

* if needed: norms depending on age and gender; also specify ethnical group, if other than caucasian; It is also accepted to refer to an already implemented analysis – however, make sure that results put together from different existing analyses are harmonized (see OnyxCeph³™ menu item Help |Findings lexicon).

For creating a customized analysis a flat rate of EUR 290 plus VAT (if applicable) is charged one-off. The analysis is delivered as a patch file to be imported via menu item Options|Install Patch. If in addition, also a print template should be created, for each page EUR 30 plus VAT (if applicable) will be charged

Further information can be found on www.onyxceph.com



Image Instruments GmbH
Niederwaldstraße 3
D 09123 Chemnitz
Germany
Tel. +49 371 9093 140
Fax +49 371 9093 149
Email info@image-instruments.de
Internet www.image-instruments.de
Register HRB 16678 (Amtsgericht Chemnitz)
Tax-No 214 111 03320
VAT-ID DE 812699756