
Required Information for Creating a Customized OnyxCeph^{3™} Analysis

The following information is required to create a customized analysis:

0. Tracing type of the analysis (Lateral Ceph, Pano, Profile Photo, Frontal Foto,)
1. Name of the analysis
2. List of all numeric results to be calculated, incl.
 - 2.1 Short name of the result variable
 - 2.2 Long name of the result variable
 - 2.3 Definition for the result variable incl. no. of digits after decimal point for the numerical value
 - 2.4 Clinical norm (if available) *
 - 2.5 Standard deviation for the norm (if available) *
 - 2.6 Verbal result description for result interbvals (if available)
3. A sketch for the graphical tracing

* if needed, clinical norms depending on age and gender; also specify ethnical group, if other than caucasian;

It is also accepted to refer to an already implemented analysis – however, make sure that results put together from different existing analyses are harmonized (see OnyxCeph^{3™} menu item Help |Findings lexicon).

For creating a customized analysis a flat rate of EUR 290 plus VAT (if applicable) is charged one-off. The analysis is delivered as a patch file to be imported via menu item Options|Install Patch. If you also order the creation of a standard template for a print form, you will be invoiced acc. To the required time. We will send you an offer in advance for confirmation.

Further information can be found on www.onyxceph.com



Image Instruments GmbH
Niederwaldstraße 3
D 09123 Chemnitz
Germany
Tel. +49 371 9093 140
Fax +49 371 9093 149
Email info@image-instruments.de
Internet www.image-instruments.de
Register HRB 16678 (Amtsgericht Chemnitz)
Tax-No 214 111 03320
VAT-ID DE 812699756